Fiche technique
User manual
Gebrauchsanweisung

17.17209
à
17.17299

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Important information for patients and their families

This cushion is a support surface designed to help preventing pressure sores. Your mobility is reduced which can result in pressures sores.

Definition of pressure sore: A pressure sore is an ischemic ulceration of the skin caused by prolonged, excess pressure of soft tissue over bony prominences of the body and a support surface. The continuous, excess pressure can hinder blood circulation and lead to the formation of a pressure sore.

There are different stages of pressure ulcer: a persistent redness despite the relieving of pressure, a blister with or without abrasion to the epidermis or dermis, a subcutaneous tissue injury that can become an ulceration extending to muscles, tendons and even bones.

Role of support surfaces: this support is designed to reduce pressure on at-risk zones by distributing body weight uniformly over the entire contact surface. Less pressure contributes to the prevention of pressure sores.

Recommendations and precautions of use: Please note that a support alone is insufficient in preventing pressure sores; it is recommended to:


- change position at least every 2 or 3 hours
- alternate with static or dynamic standing, if possible
- maintain skin hygiene and protect the skin from maceration
- in case of incontinence, change the protective elements regularly
- observe daily the condition of the skin
- ensure that patient's diet is sufficient and well-adapted
- drink regularly

If you are unable to fulfil one of these measures or in case of abnormal events, such as fever, pains, redness or whitening of weight-bearing areas (head, shoulders, back, hip, shoulder blades, pelvis, heels, ...), you have to inform the nursing staff as soon as possible.

The user of such a support as well as his/her family must be trained in the use of the device. It has to be provided by a healthcare professional at home or in a healthcare facility. A personal logbook recording recommendations and preventive actions can be created.

It is important to avoid too much thickness between the body and the support, except for underwear and possibly a complete adult brief. Use rather loose-fitting cotton underwear and if possible without seams in the weight-bearing areas if possible. Do not insert folded towel or sheet or any additional cushion between the cushion and the patient and ensure that there is no foreign bodies such as tubes, crumbs, fatty substances, ...

For hygiene reasons, each pressure sore cushion should be reserved for the exclusive use of **one patient**.

The cushion's efficiency depends on a proper positioning while sitting.


Have your seat adjusted by a professional. Relieve the pressure points as often as possible using your arms.

Important: Sliding forward promotes shearing. Incorrect adjustment of the wheelchair footrests (feet too high, too low or uneven bearing pressure), promotes pressure points.

Product description

- 1- Viscoelastic PU foam
- 2- Two Pharmatex covers
- 3- Non-slip bottom
- 4- Transport handle
- 5- Zipper

Technical characteristics

Cushion: Viscoelastic foam	Polyurethane 80 kg/m ³	
Black Pharmatex cover top	60% Polyurethane - 40% Polyester	Weight : 120 g/m ²
Non-slip cover bottom	57% Polyurethane - 43% Polyester	
Fire resistance test	EN 597-1 & EN 597-2	

This cushion meets the protocol requirements of the french regulation on bed sore prevention, and has been controlled by the FCBA which delivered conformity certificate n°ESC 19-004 from october 3, 2019.

References	Dimensions		References	Dimensions		References	Dimensions	
17.17209	34 x 36 x 8	2 covers	17.17218	40 x 42 x 8	2 covers	17.17225	45 x 40 x 8	2 covers
17.209	34 x 36 x 8	1 cover	17.218	40 x 42 x 8	1 cover	17.225	45 x 40 x 8	1 cover
17.17212	34 x 38 x 8	2 covers	17.17219	40 x 45 x 8	2 covers	17.17235	46 x 45 x 8	2 covers
17.212	34 x 38 x 8	1 cover	17.219	40 x 45 x 8	1 cover	17.235	46 x 45 x 8	1 cover
17.17214	36 x 36 x 8	2 covers	17.17228	42 x 38 x 8	2 covers	17.17229	48 x 40 x 8	2 covers
17.214	36 x 36 x 8	1 cover	17.228	42 x 38 x 8	1 cover	17.229	48 x 40 x 8	1 cover
17.17216	36 x 38 x 8	2 covers	17.17226	42 x 40 x 8	2 covers	17.17236	48 x 43 x 8	2 covers
17.216	36 x 38 x 8	1 cover	17.226	42 x 40 x 8	1 cover	17.236	48 x 43 x 8	1 cover
17.17234	37,5 x 38 x 8	2 covers	17.17295	42 x 42 x 8	2 covers	17.1727	48 x 45 x 8	2 covers
17.234	37,5 x 38 x 8	1 cover	17.295	42 x 42 x 8	1 cover	17.271	48 x 45 x 8	1 cover
17.17217	38 x 42 x 8	2 covers	17.17227	42 x 45 x 8	2 covers	17.17299	50 x 45 x 8	2 covers
17.217	38 x 42 x 8	1 cover	17.227	42 x 45 x 8	1 cover	17.299	50 x 45 x 8	1 cover
17.17294	40 x 40 x 8	2 covers	17.17267	43 x 44 x 8	2 covers	Patients up to 100 kg for widths < 45 cm 130 kg for widths > 45 cm (pressure value 169 g/cm ²)		
17.294	40 x 40 x 8	1 cover	17.267	43 x 44 x 8	1 cover			

Principles of use

The cushion is made of high-density viscoelastic foam and covered by a micro-ventilated, mite-resistant and CRIB 5 fire-resistant lined polyurethane Pharmatex cover. The non-slip bottom fabric ensures that the cushion is securely maintained on the wheelchair's seat.

The cover allows the patient to not be in direct contact with the foam for hygiene reasons and to avoid any risk of heat development.

This viscoelastic foam cushion must be used in accordance with the recommendations of a qualified healthcare professional. The level of risk for pressure sore development according to the BRADEN or NORTON scale, the weight of the user and the intended use by caregivers must be taken into consideration when selecting a suitable support.

From a technical point of view, this cushion offers the following features:

- lateral stabilization of the legs thanks to the wedge
- correct positioning of the pelvis to prevent any obliquity (the asymmetric position is a source of localized hyperpressure)
- immersion of bony prominences. This helps maintaining the continuity of blood flow in the ischial and sacro-coccygeal areas.

Installation tips

This cushion is suitable for all kinds of wheelchair whose seat sling is in a good overall condition and with no sag or deformation. The size of the cushion must be compatible with the seat to ensure a proper positioning on the wheelchair. Ensure that the top/bottom orientation of the cushion is respected. The non-slip coating must be in direct contact with the fabric or the seat on which the cushion is placed. Place the handle at the front of the cushion

Precautions for use

Only use the cover, replacement cover supplied with the cushion or protection recommended by the manufacturer to preserve the pressure sore prevention characteristics of the cushion. Do not fold the cushion during transfer operations in order to prevent any damage of the viscoelastic foam.

For user comfort, use rather cotton clothing or linens. Do not use synthetic fabric and clothing: because of frictions they can cause heats and a pelvis slip that can generate shearing on cutaneous and subcutaneous tissues.

In case of incontinence, clean the cover as required. During washing operations, fit the replacement cover on the cushion to prevent all direct contact with the cushion's surface.

Warranty

This product has a 3 years warranty from the purchase date against any manufacturing defect. The warranty doesn't include changes in colour which do not affect the product's quality. It is only valid if the product has been used accordingly to cleaning and use tips.

Cushion cleaning: saturated steam cleaning without formalin, a disinfectant product or soapy water. Allow to air dry in an open area before putting the cover back on.

Important: the protective Pharmatex cover must be replaced only once the foam is perfectly dry.

Pharmatex cover cleaning: Washing machine up to 95°C maximum, recommended up to 60°C. Dryer up to 60°C. Daily cleaning: clean with sponge and soapy water, rinse with clear water and properly dry.

Never use chlorinated products - even diluted - to clean the cover (Javel, Dakin *).

Disinfection tips: Resistance to disinfection- denatured Alcohol: 90% +2% ether, Resistance to: HAC, réocid, lyorthol
Do not use cleaning products or corrosive disinfectants such as the acetone or any industrial grease removers in order to prevent deterioration of the equipment.

Monitoring and maintenance, possible repairs, cover replacement: check the following points: holes or traces of erosion on the cover and contact your retailer in case of deterioration

Storage / Recycling

Storage: Store protected from light, heat and in a dry and flat place

Temperature conditions: MINI : -10° C / MAXI : +40° C

Degree of humidity: MINI : 30% / MAXI : 95%

Elimination of the product: do not burn or throw it in the wild, drop it at your local recycling centre.

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